



ANALIZA BIAŁEK



100003415

Laboratory of Analysis

Client: ep-company
Order number: 100003415
Name: EURO-PHARMACIES T3
Number of samples: 1 tablet

Batch number: EP-5321M4
Recive date: 2023-03-28
Start date: 2023-04-13
End date: 2023-05-10

Sample identification

Substance	Method of taking a sample	Form
Triiodothyronine	Sample provided by customer	White tablets in string bag

Results

Substance	Sample name	Content	Purity
Triiodothyronine	EURO-PHARMACIES T3 50MCG	46,42mcg/tab	-----

Research comment

The test sample was applied in the amount of 10 µl to the BioBasic 8 Thermo Scientific column. During the analysis with the method a linear gradient in the solutions of H2O and Acetonitrile with the addition of 0.1% trifluoroacetic acid, a peak of the sought substance appeared. Its concentration, when compared to the reference chromatogram, is 46,42mcg/tab. The measurement error is + - 0,5mcg. The actual state of the sample sent before the test is presented on the photo below

Photo of the test sample



Password

TWET430T

Analiza Białek sp. z o.o.
KRS: 0000966168
NIP: 8992922970
REGON: 521751875

Headquarters and Laboratory:
ul. Duńska 9, 54-427 Wrocław
Budynek DELTA

Lab technician:
Analysis and research specialist
mgr Arkadiusz Zajac




ANALIZA BIAŁEK



100003415

Chromatogram



Method

Time	H2O + 0,1% TFA	ACN + 0,1% TFA
0min	60%	40%
2min	60%	40%
27min	5%	95%
30min	5%	95%

Password

Analiza Białek sp. z o.o.
KRS: 0000966168
NIP: 8992922970
REGON: 521751875

Headquarters and Laboratory:
ul. Duńska 9, 54-427 Wrocław
Budynek DELTA

Lab technician:
Analysis and research specialist
mgr Arkadiusz Zajac