



ANALIZA BIAŁEK



100009976

## Laboratory of Analysis

**Client:** EP-COMPANY  
**Order number:** 100009976  
**Name:** EURO-PHARMACIES PRIMOBOLAN DEPOT  
**Number of samples:** 1 vial

**Batch number:** EP-7TSWBG  
**Recive date:** 2025-09-29  
**Start date:** 2025-10-10  
**End date:** 2025-10-15

## Sample identification

| Substance             | Method of taking a sample   | Form                               |
|-----------------------|-----------------------------|------------------------------------|
| Methenolone Enanthate | Sample provided by customer | Transparent liquid in a glass vial |


## Results

| Substance               | Sample name                                  | Content    | Purity |
|-------------------------|----------------------------------------------|------------|--------|
| Testosterone Propionate | EURO-PHARMACIES<br>PRIMOBOLAN DEPOT 100MG/ML | 98,57mg/ml | -----  |

## Research comment

The test sample was applied in the amount of 3µl to the BioBasic 8 Thermo Scientific column. During the course of the analysis with the method a linear gradient in the solutions of H2O and Acetonitrile with the addition of 0.1% trifluoroacetic acid, a peak of the sought substance appeared. Its concentration, when compared to the reference chromatogram, is 98,57mg/ml. Height the measurement error is + - 0,05mg. The actual state of the sample sent before the test is presented on the photo below

## Photo of the test sample



## Password

S9L5EG34

**Analiza Białek sp. z o.o.**  
**KRS:** 0000966168  
**NIP:** 8992922970  
**REGON:** 521751875

**Headquarters and Laboratory:**  
ul. Duńska 9, 54-427 Wrocław  
Budynek DELTA

**Lab technician:**  
Analysis and research specialist  
mgr Arkadiusz Zajac




ANALIZA BIAŁEK



100009976

## Chromatogram



## Method

| Time  | H2O + 0,1% TFA | ACN + 0,1% TFA |
|-------|----------------|----------------|
| 0min  | 60%            | 40%            |
| 2min  | 60%            | 40%            |
| 27min | 5%             | 95%            |
| 30min | 5%             | 95%            |

## Password

S9L5EG34

**Analiza Białek sp. z o.o.**

**KRS:** 0000966168

**NIP:** 8992922970

**REGON:** 521751875

**Headquarters and Laboratory:**

ul. Duńska 9, 54-427 Wrocław  
Budynek DELTA

**Lab technician:**

Analysis and research specialist  
mgr Arkadiusz Zajac